

Identification

Taric code: 3822 00 00

Applications

for histology, microscopy.

Specifications

appearance..... viscous, clear, colorless suitability for histology and microscopy.....passes test
liquid
refractive index at 20 °C..... 1,515 - 1,525

Physical data

- Density: 0,94 g/cm³
- Flash point: 8 °C
- Refraction index: (n 20 °C/D) 1,52

Safety - GHS

Signal Word: Danger

**Hazard Statements:**

H225: Highly flammable liquid and vapour.
H332: Harmful if inhaled.
H315: Causes skin irritation.
H361d: Suspected of damaging the unborn child.
H336: May cause drowsiness or dizziness.
H373: May cause damage to organs through prolonged or repeated exposure.
H304: May be fatal if swallowed and enters airways.

Precautionary Statements:

P210: Keep away from heat / sparks / open flames / hot surfaces. - No smoking.
P241: Use explosion-proof electrical / ventilating / lighting / equipment.
P301+P310: IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.
P303+P361+P353: IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower.
P405: Store locked up.
P501a: Dispose of contents / container in accordance with local / regional / national / international regulations.

Transport/storage

- ADR: 3 FT1 II • UN 1993 • FLAMMABLE LIQUID, N.O.S. (DPX)
- IMDG: 3 II • UN 1993 • FLAMMABLE LIQUID, N.O.S.(DPX)
- IATA/ICAO: 3 II • UN 1993 • FLAMMABLE LIQUID, N.O.S.(DPX)
- PAX: 305
- CAO: 307
- 10°C - 30°C